


## Rangap1 Rabbit Pab For Wb, Ihc-P, If/Icc, Elisa - P46060

Item Number: CM0018092



### Introduction

RanGAP1 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA. Targets human Ran GTPase-activating protein 1 (UniProt P46060), a key regulator of RAN-mediated nuclear transport. Ideal for cell biology and cancer research.

[Learn More](#)

| Parameter        | Value                                                                                                                |
|------------------|----------------------------------------------------------------------------------------------------------------------|
| Product Name     | RanGAP1 Rabbit pAb                                                                                                   |
| Remarks/Alias    | SD, Fug1, RANGAP, RanGAP1                                                                                            |
| Species          | Human                                                                                                                |
| Gene ID (Human)  | 5905                                                                                                                 |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 398-587 of human RanGAP1 (NP_002874.1) |
| Source           | Rabbit                                                                                                               |
| Category         | Polyclonal Antibodies                                                                                                |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                                             |
| Cross Reactivity | Human                                                                                                                |
| SWISS            | P46060                                                                                                               |
| Protein Weight   | 64kDa                                                                                                                |
| Shipping         | Ice bag                                                                                                              |