

Rangap1 Rabbit Pab - P46060

Item Number: CM0030807



Introduction

Rabbit polyclonal antibody against human RanGAP1 (SWISS: P46060). Validated for WB and ELISA, cross-reacts with human and mouse. Ideal for studying nucleocytoplasmic transport and GTPase activation.

[Learn More](#)

Parameter	Value
Product Name	RanGAP1 Rabbit pAb
Remarks/Alias	SD; Fug1; RANGAP; RanGAP1
Species	Human
GeneID (Human)	5905
GeneID	5905
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 398-587 of human RanGAP1 (NP_002874.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P46060
Protein Weight	64kDa
Shipping	Ice bag