


Rapgef6 Rabbit Polyclonal Antibody For Wb, Elisa - Q8Teu7

Item Number: CM0019758



Introduction

RAPGEF6 Rabbit Polyclonal Antibody validated for WB and ELISA, cross-reacts with Mouse and Rat. Ideal for studying Rap1A, Rap2A, and M-Ras GTPase signaling in brain and heart tissues. Swiss-Prot: Q8TEU7.

[Learn More](#)

Parameter	Value
Product Name	RAPGEF6 Rabbit pAb
Remarks/Alias	RAGEF2; PDZGEF2; KIA001LB; PDZ-GEF2; RA-GEF-2; RAPGEF6
Species	Human
GeneID	51735 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 145-240 of human RAPGEF6 (NP_001157858)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q8TEU7
Protein Weight	179kDa
Shipping	Ice bag