

Raptor Rabbit Polyclonal Antibody For Wb, Elisa - Q8N122

Item Number: CM0021515



Introduction

High-quality rabbit polyclonal antibody against human Raptor (mTORC1 component). Validated for WB and ELISA; reacts with human, mouse, rat. Ideal for nutrient sensing and cell growth research.

[Learn More](#)

Parameter	Value
Product Name	Raptor Rabbit pAb
Remarks/Alias	KOG1; Mip1; Raptor
Species	Human
Gene ID (human)	57521
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 888-1013 of human Raptor (NP_065812.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N122
Protein Weight	149kDa
Shipping	Ice bag