

Raptor Rabbit Polyclonal Antibody For Wb, Elisa - Q8N122

Item Number: CM0019609



Introduction

High-quality rabbit polyclonal antibody against human Raptor (RPTOR), validated for WB and ELISA. Recognizes human, mouse, and rat Raptor. Ideal for mTORC1 signaling and autophagy research.

[Learn More](#)

Parameter	Specification
Product Name	Raptor Rabbit pAb
Remarks/Alias	KOG1; Mip1; Raptor
Species	Human
GeneID	57521
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 888-1013 of human Raptor (NP_065812.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N122
Protein Weight	149kDa
Shipping	Ice bag