

Rasgrp3 Rabbit Pab - Q8lv61

Item Number: CM0026375

Introduction

RASGRP3 rabbit polyclonal antibody validated for WB, IF/ICC, and ELISA. Recognizes human and mouse RASGRP3, a calcium/DAG-regulated GEF for Ras and Rap1. Ideal for signaling studies and assay development.

[Learn More](#)

Parameter	Value
Product Name	RASGRP3 Rabbit pAb
Remarks/Alias	GRP3; RASGRP3
Species	Human
GeneID (Human)	25780
GeneID	25780
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 540-689 of human RASGRP3 (NP_056191.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8IV61
Protein Weight	78kDa
Shipping	Ice bag