


Rat Anti-Mouse Cd122 Monoclonal Antibody For Flow Cytometry - P16297

Item Number: CM0011489



Introduction

Flow cytometry-validated -conjugated rat monoclonal antibody targeting mouse CD122 (IL-2R β). Ideal for immune cell profiling and signaling studies in mouse models.

[Learn More](#)

Parameter	Value
Product Name	Rat anti-Mouse CD122 mAb
Remarks/Alias	p70; CD122; IL15Rbeta; IL-2Rbeta; IL-15Rbeta; IL-2/15Rbeta
Species	Mouse
Gene ID (Human)	3560
Gene ID (Mouse)	16185
Immunogen	Synthetic peptide derived from mouse CD122
Source	Rat
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P16297
Protein Weight	60kDa
Shipping	Ice bag