

Rat Anti-Mouse Cd3 Monoclonal Antibody For Flow Cytometry - P04235 / P22646

Item Number: CM0005688

Introduction

-conjugated rat monoclonal antibody targeting mouse CD3 delta (P04235) and epsilon (P22646) chains. Validated for flow cytometry, enabling direct detection of T cells in mouse samples. Ideal for immunophenotyping and T-cell research.

[Learn More](#)

Parameter	Specification
Product Name	Rat anti-Mouse CD3 mAb
Remarks/Alias	T3d; CD3; T3e; CD3epsilon
Species	Mouse
GeneID (Human)	916/915
GeneID	12500/12501
Immunogen	A synthesized peptide derived from mouse CD3.
Source	Rat
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P04235/P22646
Protein Weight	19kDa/21kDa
Shipping	Ice bag