

## Rat Anti-Mouse Cd45R/B220 Mab For Flow Cytometry - P06800

Item Number: CM0009147

### Introduction

Rat monoclonal antibody conjugated to targeting mouse CD45R/B220 (PTPRC, L-CA, Ly-5). Validated for flow cytometry (FC), suitable for immunophenotyping of mouse leukocyte subsets. Reliable raw material for diagnostic and research use.

[Learn More](#)

| Parameter        | Value                                      |
|------------------|--------------------------------------------|
| Product Name     | Rat anti-Mouse CD45R/B220 mAb              |
| Remarks/Alias    | B220, Cd45, L-CA, Ly-5, T200, CD45R, Lyt-4 |
| Species          | Mouse                                      |
| Gene ID (Human)  | 5788                                       |
| Gene ID          | 19264                                      |
| Immunogen        | Recombinant protein of mouse CD45R/B220    |
| Source           | Rat                                        |
| Category         | Monoclonal Antibodies                      |
| Application      | FC                                         |
| Cross Reactivity | Mouse                                      |
| SWISS            | P06800                                     |
| Protein Weight   | 130kDa/140kDa/144kDa                       |
| Shipping         | Ice bag                                    |