

Rat Anti-Mouse Cd45R/B220 Monoclonal Antibody For Fc - P06800

Item Number: CM0011705

Introduction

Rat monoclonal antibody targeting mouse CD45R/B220, suitable for flow cytometry. Recognizes the leukocyte common antigen (L-CA) expressed on hematopoietic cells.

[Learn More](#)

Parameter	Value
Product Name	Rat anti-Mouse CD45R/B220 mAb
Remarks/Alias	loc; B220; Cd45; L-CA; Ly-5; T200; CD45R; Lyt-4
Species	Mouse
Gene ID (Human)	5788
Gene ID	19264
Immunogen	Recombinant protein of mouse CD45R/B220
Source	Rat
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P06800
Protein Weight	130 kDa/140 kDa/144 kDa
Shipping	Ice bag