

Rat Anti-Mouse Cd45R/B220 Monoclonal Antibody For Flow Cytometry - P06800 / P08575

Item Number: CM0006071



Introduction

-conjugated rat monoclonal antibody specific to mouse CD45R/B220 (PTPRC, L-CA, Ly-5, T200). Validated for flow cytometry. Ideal for B-cell immunophenotyping and immune cell signaling research.

[Learn More](#)

Parameter	Value
Product Name	Rat anti-Mouse CD45R/B220 mAb
Remarks/Alias	loc; B220; Cd45; L-CA; Ly-5; T200; CD45R; Lyt-4
Species	Mouse
Gene ID (human)	5788
Gene ID	19264/5788
Immunogen	Synthetic peptide
Source	Rat
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P06800/P08575
Protein Weight	130kDa/140kDa/144kDa
Shipping	Ice bag