

Rat Anti-Mouse Cd45R/B220 Monoclonal Antibody For Flow Cytometry - P06800

Item Number: CM0005424



Introduction

Rat anti-mouse CD45R/B220 monoclonal antibody validated for flow cytometry. Targets receptor-type tyrosine-protein phosphatase C (PTPRC), a critical leukocyte marker in mouse for T-cell activation and B-cell research.

[Learn More](#)

Parameter	Value
Product Name	Rat anti-Mouse CD45R/B220 mAb
Remarks/Alias	loc; B220; Cd45; L-CA; Ly-5; T200; CD45R; Lyt-4
Species	Mouse
GenelD (Human)	5788
GenelD	19264/5788
Immunogen	Recombinant protein of mouse CD45R/B220
Source	Rat
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P06800
Protein Weight	130 kDa/140 kDa/144 kDa
Shipping	Ice bag