

Rat Anti-Mouse F4/80 Monoclonal Antibody For Fc - Q61549

Item Number: CM0007982

Introduction

High-specificity -conjugated rat monoclonal antibody targeting mouse F4/80 (EMR1), validated for flow cytometry. Ideal for macrophage identification and immune cell profiling in immunological research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	Rat anti-Mouse F4/80 mAb
Remarks/Alias	Emr1; Ly71; F4/80; Gpf480; TM7LN3; DD7A5-7; EGF-TM7; Adgre1
Species	Mouse
Gene ID (Human)	2015
Gene ID	13733
Immunogen	Recombinant protein of mouse F4/80
Source	Rat
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q61549
Protein Weight	102kDa
Shipping	Ice bag