

Rat Anti-Mouse F4/80 Monoclonal Antibody For Flow Cytometry - Q61549

Item Number: CM0012755

Introduction

Flow cytometry-validated rat anti-mouse F4/80 monoclonal antibody; reacts with mouse macrophage glycoprotein (102 kDa adhesion GPCR). Ideal for immunology studies in mouse models.

[Learn More](#)

Parameter	Value
Product Name	Rat anti-Mouse F4/80 mAb
Remarks/Alias	Emr1; Ly71; F4/80; Gpf480; TM7LN3; DD7A5-7; EGF-TM7
Species	Mouse
GenelD (human)	2015
GenelD	13733
Immunogen	Recombinant protein of mouse F4/80.
Source	Rat
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q61549
Protein Weight	102kDa
Shipping	Ice bag