

Rat Anti-Mouse F4/80 Monoclonal Antibody For Flow Cytometry - Q61549

Item Number: CM0005160

Introduction

Rat monoclonal antibody targeting mouse F4/80 (Adgre1), an adhesion GPCR on macrophages. Validated for flow cytometry (FC). This conjugate facilitates direct detection in immunological studies. Cross-reacts with mouse.

[Learn More](#)

Parameter	Value
Product Name	Rat anti-Mouse F4/80 mAb
Remarks/Alias	Emr1, Ly71, F4/80, Gpf480, TM7LN3, DD7A5-7, EGF-TM7, Adgre1
Species	Mouse
GeneID (human)	2015
GeneID	13733
Immunogen	Recombinant protein of mouse F4/80
Source	Rat
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q61549
Protein Weight	102kDa
Shipping	Ice bag