

Rat Anti-Mouse Ly6G Monoclonal Antibody For Fc - P35461

Item Number: CM0012425

Introduction

-conjugated rat monoclonal antibody against mouse Ly6G, a GPI-anchored cell membrane protein expressed in bone marrow. Designed for flow cytometry (FC) applications, enabling reliable detection in mouse samples.

[Learn More](#)

Parameter	Value
Product Name	Rat anti-Mouse Ly6G mAb
Remarks/Alias	Gr1; Gr-1; Ly-6G
Species	Mouse
Gene ID (Human)	80740
Gene ID	546644
Immunogen	Synthetic peptide A synthesized peptide derived from mouse Ly-6G.
Source	Rat
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P35461
Shipping	Ice bag