


Rat Anti-Mouse Mhc Class Ii (I-A/I-E) Monoclonal Antibody For Fc - P18468 / P06342 / Q3U060 / P14435

Item Number: CM0011300



Introduction

Rat monoclonal antibody targeting mouse MHC Class II (I-A/I-E), suitable for flow cytometry. Recognizes mouse H-2Eb, Ia4, and related alleles. Ideal for immunology research on adaptive immunity and antigen presentation.

[Learn More](#)

Parameter	Value
Product Name	Rat anti-Mouse MHC Class II (I-A/I-E) Monoclonal Antibody
Remarks/Alias	Eb; Ia4; H2Eb; Ia-4; H-2Eb
Species	Mouse
GeneID (Mouse)	14969, 14961, 381091, 14960
Immunogen	Recombinant protein of mouse MHC Class II (I-A/I-E)
Source	Rat
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P18468, P06342, Q3U060, P14435
Protein Weight	~30kDa / ~32kDa / ~26kDa
Shipping	Ice bag