


Rbfox2 Rabbit Polyclonal Antibody For Wb, Elisa - 043251

Item Number: CM0021431



Introduction

High-quality rabbit polyclonal antibody against RBFOX2 (RNA binding protein fox-1 homolog 2), validated for Western blot (WB) and ELISA with human reactivity. Ideal for alternative splicing, mRNA processing, and gene regulation studies.

[Learn More](#)

Parameter	Value
Product Name	RBFOX2 Rabbit pAb
Remarks/Alias	RTA; fxh; FOX2; RBM9; Fox-2; HNRBP2; HRNBP2; dj106l20.3; RBFOX2
Species	Human
GeneID (Human)	23543
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-110 of human RBFOX2 (NP_001026865.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O43251
Protein Weight	41kDa
Shipping	Ice bag