

Rbm20 Rabbit Polyclonal Antibody For Wb, Elisa - Q5T481

Item Number: CM0024330

Introduction

RBM20 Rabbit Polyclonal Antibody validated for WB and ELISA, targeting human RBM20, a cardiac splicing regulator. Ideal for cardiovascular disease research. - Q5T481

[Learn More](#)

Parameter	Value
Product Name	RBM20 Rabbit pAb
Species	Human
Gene ID	282996 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 820-980 of human RBM20 (NP_001127835.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q5T481
Protein Weight	134kDa
Shipping	Ice bag