

Rbmx/Hnrnp G Rabbit Mab For Wb, Ihc-P, Elisa - P38159

Item Number: CM0003243



Introduction

Rabbit monoclonal antibody against human RBMX/hnRNP G (SWISS P38159), suitable for WB, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat. Targets a key RNA-binding protein involved in alternative splicing and tumor suppression.

[Learn More](#)

Parameter	Value
Product Name	RBMX/hnRNP G Rabbit mAb
Remarks/Alias	RNMx; HNRPG; HNRNPG; MRXS11; RBMXP1; RBMXRT; hnRNP-G; RBMX/hnRNP G
Species	Human
GeneID (Human)	27316
GeneID	27316
Immunogen	A synthetic peptide corresponding to amino acids 292-391 of human RBMX/hnRNP G (NP_002130.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P38159
Protein Weight	42kDa
Shipping	Ice bag