


Rbmy1A1 Rabbit Pab - P0Djd3

Item Number: CM0019373



Introduction

RBMY1A1 Rabbit Polyclonal Antibody validated for IF-P, IHC-P, ELISA. Reacts with human, mouse, rat. Targets 56kDa testis-specific splicing factor RBMY1A1 (SWISS: P0DJd3).

[Learn More](#)

Parameter	Value
Product Name	RBMY1A1 Rabbit pAb
Remarks/Alias	RBM; RBM1; RBM2; RBMY; YRRM1; YRRM2; RBMY1C; RBMY1A1
Species	Human
Gene ID (Human)	5940
Gene ID	5940
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 90-300 of human RBMY1A1 (NP_005049.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P0DJd3
Protein Weight	56kDa
Shipping	Ice bag