


Rbp2 Rabbit Pab - P50120

Item Number: CM0024749



Introduction

CamelBio's RBP2 Rabbit pAb is a rabbit polyclonal antibody against human RBP2 (CRBP-II), validated for WB, IF/ICC, ELISA in human, mouse, rat. Facilitates intracellular retinol transport studies.

[Learn More](#)

Parameter	Value
Product Name	RBP2 Rabbit pAb
Remarks/Alias	CRBP2; RBPC2; CRBP2; CRABP-II; RBP2
Species	Human
GeneID (human)	5948
GeneID	5948
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-134 of human RBP2 (NP_004155.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50120
Protein Weight	16kDa
Shipping	Ice bag