


Rbx1 Rabbit Pab For Wb And Elisa - P62877

Item Number: CM0017939



Introduction

Rabbit polyclonal antibody targeting human RBX1 (ROC1), validated for WB and ELISA, cross-reacts with mouse and rat. Ideal for studying ubiquitin ligase function in cell cycle and DNA repair pathways.

[Learn More](#)

Parameter	Value
Product Name	RBX1 Rabbit pAb
Remarks/Alias	ROC1; RNF75; BA554C12.1; RBX1
Species	Human
GeneID (Human)	9978
GeneID	9978
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human RBX1 (NP_055063.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P62877
Protein Weight	12kDa
Shipping	Ice bag