


Rchy1 Rabbit Monoclonal Antibody For Wb And Elisa - Q96Pm5

Item Number: CM0003120



Introduction

RCHY1 rabbit monoclonal antibody for WB and ELISA, targeting human RCHY1 (Q96PM5), an E3 ubiquitin ligase involved in p53 degradation and ribosomal quality control. Ideal for cancer and ubiquitination research.

[Learn More](#)

Parameter	Details
Product Name	RCHY1 Rabbit mAb
Remarks/Alias	ZCHY; ARNIP; CHIMP; PIRH2; RNF199; ZNF363; PRO1996; RCHY1
Species	Human
Gene ID	25898 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 160-261 of human RCHY1 (Q96PM5)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q96PM5
Protein Weight	30kDa
Shipping	Ice bag