

Rec8 Rabbit Polyclonal Antibody (Mouse) - Rec8

Item Number: CM0016701

Introduction

Rabbit polyclonal antibody targeting mouse REC8 protein (67kDa). Validated for Western blot (WB) and ELISA applications, with cross-reactivity to rat. Recombinant immunogen corresponds to amino acids 194-561. Ideal for meiosis and chromosome cohesion research.

[Learn More](#)

Parameter	Value
Product Name	REC8 Rabbit pAb
Remarks/Alias	mrec; Rec8L1; REC8
Species	Mouse
Gene ID (Human)	9985
Gene ID (Mouse)	56739
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 194-561 of mouse REC8 (NP_064386.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
Protein Weight	67kDa
Shipping	Ice bag