

Reg3G Rabbit Pab - Mouse Reg3G

Item Number: CM0019969

Introduction

Polyclonal antibody targeting mouse REG3G (REG-3-gamma) with specificity for human, mouse, and rat samples. Suitable for Western blot, immunofluorescence/immunocytochemistry, and ELISA applications in tissue repair and inflammatory disease research.

[Learn More](#)

Parameter	Value
Product Name	REG3G Rabbit pAb
Remarks/Alias	REG-3-gamma; reg III-gamma; REG3G
Species	Mouse
Gene ID (Human)	130120
Gene ID	19695
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 27-174 of mouse REG3G (NP_035390.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
Protein Weight	19kDa
Shipping	Ice bag