

Rel Rabbit Pab For Wb, Elisa - Q04864

Item Number: CM0017617

Introduction

Rabbit polyclonal antibody targeting human REL (c-Rel) for Western blot (WB) and ELISA applications. Reacts with human, mouse, and rat samples. Involved in NF-kB signaling and immune regulation.

[Learn More](#)

Parameter	Value
Product Name	REL Rabbit pAb
Remarks/Alias	C-Rel; IMD92; HIVEN86A; REL
Species	Human
GeneID (Human)	5966
GeneID	5966
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human REL (NP_002899.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q04864
Protein Weight	69kDa
Shipping	Ice bag