

Rev-Erb α /Nr1D1 Rabbit Pab - P20393

Item Number: CM0023846



Introduction

Rev-Erb α /NR1D1 Rabbit Polyclonal Antibody validated for WB, IHC-P, IF/ICC, ELISA in human, mouse, rat. Explore circadian rhythm, metabolism, and inflammation pathways.

[Learn More](#)

Parameter	Value
Product Name	Rev-Erb α /NR1D1 Rabbit pAb
Remarks/Alias	EAR1; hRev; THRA1; THRAL; ear-1; REVERBA; REVERBalpha; Rev-Erb α /NR1D1
Species	Human
GeneID (human)	9572
GeneID	9572
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 280-360 of human Rev-Erb α /NR1D1 (NP_068370.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20393
Protein Weight	67kDa
Shipping	Ice bag