

Rfc2 Rabbit Polyclonal Antibody For Wb And Elisa - P35250

Item Number: CM0026696



Introduction

High-quality rabbit polyclonal antibody against human RFC2 (Replication factor C subunit 2). Validated for WB and ELISA. Cross-reacts with human and mouse RFC2. Ideal for DNA replication and PCNA loading studies.

[Learn More](#)

Parameter	Value
Product Name	RFC2 Rabbit pAb
Remarks/Alias	RFC40; RFC2
Species	Human
GeneID (human)	5982
GeneID	5982
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 38-297 of human RFC2 (NP_852136.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P35250
Protein Weight	39kDa
Shipping	Ice bag