

Rhof Rabbit Pab For Wb And Elisa - Q9Hbh0

Item Number: CM0027150

Introduction

Rabbit polyclonal antibody against human RHOF (Rif), validated for WB and ELISA. Cross-reacts with human and mouse. Ideal for studying filopodia formation and actin cytoskeleton dynamics.

[Learn More](#)

Parameter	Value
Product Name	RHOF Rabbit pAb
Remarks/Alias	RIF; ARHF; RHOF
Species	Human
Gene ID (Human)	54509
Gene ID	54509
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-120 of human RHOF (NP_061907.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9HBH0
Protein Weight	24kDa
Shipping	Ice bag