


Rhou Rabbit Pab - Q7L0Q8

Item Number: CM0018900



Introduction

RHO Rabbit pAb is a polyclonal antibody validated for WB and ELISA, targeting human RHO with cross-reactivity with mouse. Ideal for cell migration and cytoskeleton research.

[Learn More](#)

Parameter	Value
Product Name	RHO Rabbit pAb
Remarks/Alias	ARHU; G28K; WRCH1; hG28K; CDC42L1; RHO
Species	Human
GeneID (Human)	58480
GeneID	58480
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 40-100 of human RHO (NP_067028.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q7L0Q8
Protein Weight	28kDa
Shipping	Ice bag