

Rig-I/Ddx58 Rabbit Mab For Wb, Elisa - O95786

Item Number: CM0005358



Introduction

Monoclonal rabbit anti-Rig-I/DDX58 antibody validated for Western blot and ELISA. Detects human and mouse RIG-I, an innate immune receptor for viral RNA, with recombinant protein immunogen.

[Learn More](#)

Parameter	Value
Product Name	Rig-I/DDX58 Rabbit mAb
Remarks/Alias	RIG1; DDX58; RIG-I; RLR-1; SGMRT2; Rig-I/DDX58
Species	Human
GeneID (Human)	23586
GeneID	23586
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O95786
Protein Weight	107kDa
Shipping	Ice bag