

Ring1A Polyclonal Antibody (Rabbit) For Wb, Ihc-P, Elisa - Q06587

Item Number: CM0022393

Introduction

Rabbit polyclonal antibody to RING1A (RNF1, RING1) for Western blot, IHC-P, and ELISA. Validated in human, mouse, and rat. Targets E3 ubiquitin-protein ligase RING1 involved in histone H2A monoubiquitination and transcriptional repression.

[Learn More](#)

Parameter	Value
Product Name	RING1A Rabbit pAb
Remarks/Alias	RNF1; RING1A
Species	Human
Gene ID (Human)	6015
Gene ID	6015
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 207-406 of human RING1A A (NP_002922.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q06587
Protein Weight	42kDa
Shipping	Ice bag