

Ring1A Rabbit Mab For Wb, Ihc-P, Elisa - Q06587

Item Number: CM0002294

Introduction

High-quality rabbit mAb against human RING1A, validated for WB, IHC-P, ELISA. Targets E3 ligase involved in H2A ubiquitination and transcriptional repression. Ideal for epigenetics and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	RING1A Rabbit mAb
Remarks/Alias	RNF1; RING1A
Species	Human
GenelD (Human)	6015
GenelD	6015
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q06587
Protein Weight	42kDa
Shipping	Ice bag