

Rnf10 Rabbit Polyclonal Antibody For Wb, Elisa - Q8N5U6

Item Number: CM0017769

Introduction

Rabbit polyclonal antibody against human RNF10, raised against a recombinant fusion protein (aa 552-811). Validated for WB and ELISA, cross-reacts with mouse and rat. Suitable for studying ribosome quality control and ubiquitination pathways.

[Learn More](#)

Parameter	Value
Product Name	RNF10 Rabbit pAb
Remarks/Alias	RIE2; RNF10
Species	Human
GeneID (Human)	9921
GeneID	9921
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q8N5U6
Protein Weight	90kDa
Shipping	Ice bag