

Rnf138 Rabbit Pab - Rnf138

Item Number: CM0018692

Introduction

Rabbit polyclonal antibody targeting RNF138 for WB and ELISA, with verified cross-reactivity to human. Suitable for ubiquitination and protein degradation research.

[Learn More](#)

Parameter	Value
Product Name	RNF138 Rabbit pAb
Remarks/Alias	Trif; STRIN; Trif-d; 2410015A17Rik; 2810480D20Rik; RNF138
Species	Mouse
Gene ID (Human)	51444
Gene ID (Mouse)	56515
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 60-245 of mouse RNF138 (Q9CQE0)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
Protein Weight	24kDa/28kDa
Shipping	Ice bag