


Rnf20 Rabbit Polyclonal Antibody For Wb, Elisa - Q5Vtr2

Item Number: CM0022957



Introduction

High-quality rabbit polyclonal antibody against human RNF20 (BRE1A), validated for Western blot and ELISA. Reacts with human, mouse, rat. Suitable for studying histone H2B ubiquitination and chromatin regulation. - Q5VTR2

[Learn More](#)

Parameter	Value
Product Name	RNF20 Rabbit pAb
Remarks/Alias	BRE1; BRE1A; hBRE1; RNF20
Species	Human
GeneID (Human)	56254
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human RNF20 (NP_062538.5).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5VTR2
Protein Weight	114kDa
Shipping	Ice bag