

## Romo1 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P60602

Item Number: CM0019512

### Introduction

Reactive oxygen species modulator 1 (ROMO1) rabbit polyclonal antibody validated for WB, IF-P, ELISA; cross-reacts with human, mouse, rat; mitochondrial inner membrane protein involved in ROS production and cancer.

[Learn More](#)

| Parameter        | Value                                                                                                                  |
|------------------|------------------------------------------------------------------------------------------------------------------------|
| Product Name     | ROMO1 Rabbit pAb                                                                                                       |
| Remarks/Alias    | MTGM; MTGMP; C20orf52; bA353C18.2; ROMO1                                                                               |
| Species          | Human                                                                                                                  |
| GeneID (human)   | 140823                                                                                                                 |
| GeneID           | 140823                                                                                                                 |
| Immunogen        | Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-79 of human ROMO1(NP_542786.1). |
| Source           | Rabbit                                                                                                                 |
| Category         | Polyclonal Antibodies                                                                                                  |
| Application      | WB, IF-P, ELISA                                                                                                        |
| Cross Reactivity | Human,Mouse,Rat                                                                                                        |
| SWISS            | P60602                                                                                                                 |
| Protein Weight   | 8kDa                                                                                                                   |
| Shipping         | Ice bag                                                                                                                |