


Ror1 Rabbit Mab For Fc - Q01973

Item Number: CM0006479



Introduction

ROR1 Rabbit monoclonal antibody for flow cytometry (FC), targeting human ROR1 (104 kDa). Recombinant immunogen, cross-reactive with human only. Suitable for oncology and auditory research.

[Learn More](#)

Parameter	Value
Product Name	ROR1 Rabbit mAb
Remarks/Alias	NTRKR1; dj537F10.1; ROR1
Species	Human
GeneID (Human)	4919
GeneID	4919
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-403 of Human ROR1 (NP_005003.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q01973
Protein Weight	104kDa
Shipping	Ice bag