


Rpap3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9H6T3

Item Number: CM0015084



Introduction

Rabbit polyclonal antibody against human RPAP3, validated for WB, IF/ICC, and ELISA in human, mouse, and rat. RPAP3 links RNA polymerase II to chaperones for protein complex assembly.

[Learn More](#)

Parameter	Specification
Product Name	RPAP3 Rabbit pAb
Remarks / Alias	Tah1; hSpagh; RPAP3
Species	Human
Gene ID	79657
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 17-241 of human RPAP3 (NP_001139547)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H6T3
Protein Weight	76kDa
Shipping	Ice bag