

Rpb3/Polr2C Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P19387

Item Number: CM0010027



Introduction

Rabbit monoclonal antibody targeting human RPB3/POLR2C (P19387), validated for WB, IHC-P, ELISA, with cross-reactivity to mouse and rat. Ideal for studying RNA polymerase II transcription machinery.

[Learn More](#)

Parameter	Value
Product Name	RPB3/POLR2C Rabbit mAb
Remarks/Alias	RPB3; RPB31; hRPB33; hsRPB3; RPB3/POLR2C
Species	Human
Gene ID (Human)	5432
Gene ID	5432
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P19387
Protein Weight	31kDa
Shipping	Ice bag