


Rpl28 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P46779

Item Number: CM0015550



Introduction

Rabbit polyclonal antibody targeting human RPL28 (60S ribosomal protein L28). Validated for Western blot, IF/ICC, and ELISA. Cross-reacts with human, mouse, rat. Molecular weight ~16 kDa. Suitable for ribosome research.

[Learn More](#)

Parameter	Value
Product Name	RPL28 Rabbit pAb
Remarks/Alias	L28; eL28; RPL28
Species	Human
GeneID (Human)	6158
GeneID	6158
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human RPL28 (NP_000982.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P46779
Protein Weight	16kDa
Shipping	Ice bag