


Rpl35A Rabbit Polyclonal Antibody For Wb And Elisa - P18077

Item Number: CM0023855



Introduction

Rabbit polyclonal antibody targeting human RPL35A, a large ribosomal subunit protein eL33, validated for WB and ELISA with cross-reactivity to mouse and rat. Supports ribosome and hematopoietic research.

[Learn More](#)

Parameter	Value
Product Name	RPL35A Rabbit pAb
Remarks/Alias	DBA5; L35A; eL33; RPL35A
Species	Human
GeneID (Human)	6165
GeneID	6165
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 10-90 of human RPL35A (NP_000987.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P18077
Protein Weight	13kDa
Shipping	Ice bag