

Rpl7 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P18124

Item Number: CM0027478



Introduction

RPL7 Rabbit Polyclonal Antibody validated for WB, IHC-P, IF/ICC, ELISA. Recognizes human, mouse, rat RPL7. Research-use tool for ribosomal protein studies and translation apparatus analysis.

[Learn More](#)

Parameter	Value
Product Name	RPL7 Rabbit pAb
Remarks/Alias	L7; uL30; humL7-1; RPL7
Species	Human
Gene ID (Human)	6129
Gene ID	6129
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-248 of human RPL7 (NP_000962.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P18124
Protein Weight	29kDa
Shipping	Ice bag