


Rps15A Rabbit Pab For Wb, Ihc-P, Elisa - P62244

Item Number: CM0017241



Introduction

Rabbit polyclonal antibody targeting human RPS15A (uS8), suitable for WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Useful for studying small ribosomal subunit function, erythropoiesis, and Diamond-Blackfan anemia.

[Learn More](#)

Parameter	Value
Product Name	RPS15A Rabbit pAb
Remarks/Alias	uS8; S15a; DBA20; RPS15A
Species	Human
GeneID (human)	6210
GeneID	6210
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human RPS15A (NP_001010.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62244
Protein Weight	15kDa
Shipping	Ice bag