


Rps20 Rabbit Pab For Wb, If/Icc, Elisa - P60866

Item Number: CM0013497



Introduction

RPS20 Rabbit pAb is a polyclonal antibody validated in WB, IF/ICC, and ELISA for detecting RPS20 in human, mouse, and rat samples. A reliable tool for ribosomal protein research.

[Learn More](#)

Specification	Details
Product Name	RPS20 Rabbit pAb
Remarks/Alias	S20; uS10; RPS20
Species	Human
GeneID (Human)	6224
GeneID	6224
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-119 of human RPS20 (NP_001014.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P60866
Protein Weight	13kDa
Shipping	Ice bag