


Rps3A Rabbit Pab For Wb, Ihc-P, If/Icc, Ip, Elisa - P61247

Item Number: CM0025474



Introduction

Rabbit polyclonal antibody targeting RPS3A (40S ribosomal protein S3a). Suitable for WB, IHC-P, IF/ICC, IP, ELISA in human, mouse, rat. Ideal for ribosome biogenesis and protein synthesis studies.

[Learn More](#)

Parameter	Value
Product Name	RPS3A Rabbit pAb
Remarks/Alias	S3A; eS1; FTE1; MFTL; RPS3A
Species	Human
GeneID (Human)	6189
GeneID	6189
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-264 of human RPS3A (NP_000997.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61247
Protein Weight	30kDa
Shipping	Ice bag