


Rrp9 Rabbit Pab For Wb, Elisa - O43818

Item Number: CM0019146



Introduction

RRP9 Rabbit Polyclonal Antibody (O43818) for Western blot and ELISA applications. Detects human and mouse RRP9, a U3 snoRNP-associated protein involved in pre-rRNA processing.

[Learn More](#)

Parameter	Value
Product Name	RRP9 Rabbit pAb
Remarks/Alias	U3-55K; RNU3IP2; RRP9
Species	Human
GeneID (Human)	9136
GeneID	9136
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 156-475 of human RRP9 (NP_004695.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	O43818
Protein Weight	52kDa
Shipping	Ice bag