

Runx1 Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - Q01196

Item Number: CM0022283



Introduction

RUNX1 rabbit polyclonal antibody validated for WB, IF/ICC, IHC-P, ELISA. Reacts with human, mouse, rat RUNX1 transcription factor essential for hematopoiesis and leukemia research.

[Learn More](#)

Parameter	Value
Product Name	RUNX1 Rabbit pAb
Remarks/Alias	AML1; CBFA2; EVI-1; AMLCR1; PEBP2aB; CBF2alpha; AML1-EVI-1; PEBP2alpha; RUNX1
Species	Human
GeneID (Human)	861
Immunogen	Synthetic peptide (amino acids 171-260 of human RUNX1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q01196
Protein Weight	49 kDa
Shipping	Ice bag