

Runx1/2/3 Rabbit Monoclonal Antibody For Wb, Elisa - Q01196 / Q13761 / Q13950

Item Number: CM0009100



Introduction

RUNX1/2/3 rabbit monoclonal antibody for Western blot (WB) and ELISA, recognizing human and mouse RUNX1/2/3. Essential for hematopoiesis and T-cell development. Ideal for cancer research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	RUNX1/2/3 Rabbit mAb
Remarks/Alias	AML1; AML1-EVI-1; AMLCR1; CBF2alpha; CBFA2; EVI-1; PEBP2aB; PEBP2alpha; RUNX1/2/3
Species	Human
Gene ID (Human)	861/864/860
Gene ID	861/864/860
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 354-453 of human RUNX1/2/3 (NP_001001890.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q01196/Q13761/Q13950
Protein Weight	51kDa
Shipping	Ice bag